



Make Every Minute Count

AFINION™ 2 ANALYSER

The Afinion™ 2 Analyser is the result of years of continuous product innovation and refinement and delivers actionable point-of-care test results in just minutes.



www.ethitech.co.za

021 855 0307 | info@ethitech.co.za

Distributed by:

ethitech
Health Innovation

AFINION™ 2 ANALYSER

With its compact size and panel of tests, the Afinion™ 2 System is ideal for point-of-care testing in physician offices, clinics, community health centers, retirement homes, emergency rooms and hospital outpatient clinics. From just a fingerstick-sized sample, highly accurate results for CRP, HbA1c, ACR and Lipid Panel are made available during the consultation.

BENEFITS TO POINT-OF-CARE TESTING

Reduced workload in the practice¹:

- Fewer follow-up calls
- Fewer follow-up letters
- Fewer patient revisits

More informed decision-making improves patient-clinician relationships and patient motivation and attitude.²

TESTS

Afinion™ tests can be used with both the Afinion™ 2 Analyser and the Afinion™ AS100 Analyser.

PRODUCT NAME	PRODUCT CODE
AFINION™ 2 ANALYSER	1116770
AFINION™ HbA1c	1116795
AFINION™ ACR	1116781
AFINION™ CRP	1116787
AFINION™ LIPID PANEL	1116801

TECHNICAL INFORMATION

- Size: 200 mm W × 186 mm H × 328 mm D
- Weight: 3.4 kg
- LCD colour touch display
- Patient and QC ID and results stored in memory
- No maintenance required
- No operator calibration
- Extensive and comprehensive fail-safe systems prevent the possibility of erroneous results
- QC and operator lockout

CONNECTIVITY

The Afinion 2 Analyser can reliably transfer test results to an information system. Use the Ethernet® cable to interface the analyser to an information system. The Afinion 2 Analyser connectivity automatically transfers patient and control results to a connected LIS/HIS/EMR system via TCP/IP networking using the protocols POCT1-A, HL7®, ASTM

1381-95 (low level) or ASTM 1394-97 (high level), selectable by configuration. With built-in connectivity, the Afinion 2 Analyser can communicate with LIS and HIS systems to share and manage data more effectively, contributing to better patient management everywhere.

3-STEP PROCEDURE

- 1 Collect the sample with the integrated sampling device.



- 2 Place the sampling device back in the test cartridge.



- 3 Place the test cartridge in the analyser and close the lid. The processing starts automatically.



1. Crocker JB, Lee-Lewandrowski E, Lewandrowski N, et al. Implementation of point-of-care testing in an ambulatory practice of an academic medical center. *Am J Clin Pathol.* 2014;142:640-6.

2. Laurence CO, Gialamas A, Bubner T. Patient satisfaction with point-of-care testing in general practice. *Br J Gen Pract.* 2010;60(572):e98-e104.

© 2020 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. Any photos displayed are for illustrative purposes only. Any person depicted in such photos is a model. 10003944-04 10/20